



- No Messy Impressions
- Save in labor and material costs
- Decrease chair time for patients
- Better fit and improvement of comfort
- Avoid additional doctor visit when ordering additional retainers
- Give your patients the peace of mind on potentially lost or damaged retainers

clique Lite	clique Plus	clique Pro
<ul style="list-style-type: none"> • 4 Essix Plus Retainers at Debond • 3D Printed Model After Debond • Digitally Stored Orthodontic Records • Ability to remake retainers immediately if needed • Brackets Can be Digitally Removed 	<ul style="list-style-type: none"> • 4 Essix Plus Retainers at Debond • 4 Additional Essix Plus Retainers as needed over 3 year period • 3D Printed Model Before Bonding • 3D Printed Model After Debond • Digitally Stored Orthodontic Records • Ability to remake retainers immediately if needed • Brackets Can be Digitally Removed 	<ul style="list-style-type: none"> • Indirect Bonding trays • 4 Essix Plus Retainers at Debond • 4 Additional Essix Plus Retainers as needed over 3 year period • 3D Printed Model Before Bonding • 3D Printed Model After Debond • Digitally Stored Orthodontic Records • Ability to remake retainers immediately if needed • Brackets Can be Digitally Removed

